# *Avicenum*

## **ANTI-TROMBO**

### **PREMIUM**

#### **(GB) Instructions for Use**

Avicenum ANTI-TROMBO PREMIUM Calf-length stockings  
Avicenum ANTI-TROMBO PREMIUM Thigh-high stockings

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Avicenum ANTI-TROMBO PREMIUM Wadenstrümpfe  
Avicenum ANTI-TROMBO PREMIUM Schenkelstrümpfe

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# (GB)

## INSTRUCTION MANUAL

Dear customer, In your hands you are now holding compression stockings that will help you in the treatment and prevention of venous or lymphatic diseases. You can maximize the benefits of the product by its proper use and maintenance.

**Avicenum ANTI-TROMBO PREMIUM** are compression stockings of compression class I. (CCL 1). The product was manufactured by the circular knitting technology. The materials used for the manufacture of the product are non-irritating and latex-free. The mechanism behind the effects of the product is based on the action of precisely defined graduated pressure on the leg, which reduces from the ankle up to the heart. The pressure applied in the ankle area in point b (see the size table) is the greatest and reaches 18–21 mmHg. In order to achieve the desired effects, you must select the correct size of the product. These compression stockings are primarily intended for bed-bound patients and serve as a medical device for the prevention of the thromboembolic disease. There is an inspection hole located under the toes which allows checking of the correct compression efficiency (by checking for cyanosis or cold toes). **Avicenum ANTI-TROMBO PREMIUM** serve as an alternative to classic bandages and are a medical device intended for use both at home and in medical facilities.

### Product name:

## **Avicenum ANTI-TROMBO PREMIUM Calf-length stockings**

### Product variants:

- inspection hole, CCL I, sizes S – XXL

## **Avicenum ANTI-TROMBO PREMIUM Thigh-high stockings**

### Product variants:

- hem, inspection hole, CCL I, sizes S – XXL

### User profile:

**Avicenum ANTI-TROMBO PREMIUM** is a medical device intended for use both at home and in medical facilities.

### Intended use:

The product is intended for the compression therapy of venous conditions and swelling of the legs of various origin.

### Indications:

- chronic venous insufficiency – subjective and objective manifestations, prevention
- primary and secondary varices of the legs
- superficial venous thrombophlebitis
- swelling associated with post-thrombotic syndrome after deep venous phlebothrombosis
- swelling and varices during pregnancy
- conditions after phlebothrombosis and thrombophlebitis of the legs
- prevention of deep venous phlebothrombosis in surgery and associated thromboembolic diseases (especially surgery, orthopaedics, gynaecology)
- swelling of the legs (post-operative and post-traumatic, lymphatic, overweight swelling, idiopathic swelling of the legs)

- compression support after sclerotherapy of leg varices
- compression support after surgical extirpation of varices (standard surgery, radiofrequency ablation, or laser ablation)
- travel thrombosis (economy class syndrome)

If this is the first time you use compression medical hosiery, consult your doctor to select the proper compression class of the product or have your medical condition regularly examined.

#### **Contraindications:**

- acute weeping skin manifestations
- acute limb ischemia
- critical limb ischemia
- decompensated cardiac insufficiency
- phlegmasia coerulea dolens, phlegmasia alba dolens
- hypersensitivity to materials used
- limb sensitivity disorders (e.g., peripheral neuropathy in diabetes mellitus)

Your doctor should consider the use of the product in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

#### **Target group:**

Patients with primary varices, patients after venous thrombosis, and patients with the above indications. **Avicenum ANTI-TROMBO PREMIUM** prevent the development of the thromboembolic disease after surgery. The product may be used only on healthy and intact skin. The product is intended for all age groups and both sexes. The product may be used in pregnancy without limitation.

#### **Side effects:**

Side effects may occur especially when using a wrong size of the product. This may cause pressure sores, impaired blood circulation due to compression of blood vessels, or deterioration of sensitivity due to excessive compression of the nerves. Sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) at the place of contact of the product materials with the skin.

#### **How to choose the correct size:**

In order to achieve the desired effects, you must select the correct size of the product. Determine the size according to the size table by measuring the circumference at the points indicated in the figure.

#### **Note:**

Any leg swelling should be examined by your doctor. If this is the first time you experience the swelling, consult your doctor before using the product. Do not use the product if you have not determined the cause of the swelling. Always consult your doctor before combining the product with other medical devices (e.g., sleeves or braces).


- Use a tape measure.
- Take the measurements standing up, preferably in the morning right after waking up – the legs may swell during the day.

- Get assistance from another person.
- Measured values may change during the therapy – always verify the correct size before using a new product.

#### SIZE TABLE:

SIZE	Circumference in measured point		
	<b>b</b>	<b>d</b>	<b>g</b>
<b>S</b>	17 – 20	28 – 34	41 – 54
<b>M</b>	20 – 23	31 – 37	48 – 61
<b>L</b>	23 – 26	34 – 40	55 – 68
<b>XL</b>	26 – 29	37 – 43	62 – 75
<b>XXL</b>	29 – 32	40 – 46	69 – 82

- b** circumference above the ankle
- d** circumference below the knee
- g** thigh circumference (5 cm under the crotch)

If the measured values do not correspond to any of the sizes in the table (e.g., point b corresponds to size S but point d is already size L), visit a specialized medical supplies store or the manufacturer.

If the measured values differ significantly for the right and the left leg, visit a specialized medical supplies store or directly the manufacturer.

If you are on the border between two sizes, choose the bigger one. If you have difficulties choosing the size of the product, visit a specialized medical supplies store or directly the manufacturer.

If you have difficulties putting on the product, visit a specialized medical supplies store and request a demonstration of the correct method of applying it. You can also use one of the aids for putting on the product.

If the hem does not hold on the leg (in thigh-high stockings) clean the silicone straps or studs with an alcohol solution or shave the hair at the point of contact of the skin with the silicone straps or studs.

#### **With a correctly selected size:**

- you should feel noticeable but not unpleasant pressure over the entire length of the product
- the product does not slide nor roll down below the knee or on the instep when walking
- the hem does not choke (in thigh-high stockings)

Information on the compression class, size, and manufacturer can be found on the tag sewn at the top edge of product. If the tag is removed over the course of using the product, its size is still indicated on the sole of the product.

### Method of use:

Put on the product by gradually pulling it over the leg (do not pull only the top edge). It is preferable to put on the product in the morning before any physical strain or formation of swelling. When putting on, the product may not come into contact with sharp objects (e.g., jewellery) or long nails. Intensive care for your feet (no rough skin on the heels) will decrease the risk of damaging the product and increase its life. Stretch the product evenly without creases for effective distribution of pressure. If you use body cosmetics, wait about 15 minutes after application before you put on the product. In the case of visible mechanical damage to the product, stop wearing it and preferably replace it with a new product.

### Care and maintenance:

**Avicenum ANTI-TROMBO PREMIUM Thigh-High Stockings** are intended for repeated use by a single person. The product will maintain its properties and safety for six months from the first use under these conditions:

- the product can be washed in a washing machine (maximum washing temperature 95 °C, gentle cycle) in a soapy solution or using a special preparation for washing compression hosiery using a protective container (wash bag)
- do not use fabric softeners
- squeeze the wet product between two towels and dry in a horizontal position away from a direct heat source (radiator, sun) or you can dry the product in a tumble dryer at low temperature (at no more than 60 °C)
- do not use centrifugation, squeeze the wet product between two towels
- dry the product in a horizontal position away from a direct heat source (radiator, sun)
- do not iron the product
- the maximum number of washing cycles is 120
- store in a dry and dark place, preferably in the original packaging, and do not expose to direct sunlight
- do not interfere mechanically in the product, do not repair places damaged by use, and do not cut the hem or lace of the product (in thigh-high stockings)
- you can restore the adhesive properties of the silicon studs of the hem (in thigh-high stockings) by applying an alcohol solution
- the product may not come into contact with organic solvents or bleaching agents

The product is safe to use for 60 months from the date of manufacture. This date is included in the product's batch number in the format YY/MM-1234567 where YY is the year and MM the month of manufacture.

### Care symbols:

*Do not use fabric softeners*

### Prohibited use:

Stop using the product immediately if:

- if it shows signs of mechanical damage (replace the product)
- your health has deteriorated which could be associated with the use of the product – consult your doctor

**Composition:**

Calf-length stockings: 70 % PAD Nylon, 30 % elastan LYCRA®

Thigh-high stockings: 70 % PAD Nylon, 30 % elastan LYCRA®, silicon

This medical device does not contain an active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, tissues or cells of animal origin or derivatives thereof.

**Disposal:**

The product can be disposed of with normal household waste or in a textile waste container. Although the product is intended for use on healthy and intact skin, it may be contaminated with body fluids during use. In that case, dispose of the product as an infectious material.

**The material used is not toxic unless burned. Never dispose of the product by burning it. Any serious adverse event that occurred in connection with the use of this product should be reported to the manufacturer and the competent public authorities of the user's country of residence.**

Czech Republic	Státní ústav pro kontrolu léčiv, Šrobárova 48, 100 41 Prague 10
Slovak Republic	Štátny ústav pre kontrolu liečiv, Kvetná 1024/11, 821 08 Bratislava
Poland	Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych, Al. Jerozolimskie 181C, 02-222 Warsaw
Germany	Das Bundesinstitut für Arzneimittel und Medizinprodukte, Kurt Georg Kiesinger Allee 3, 53175 Bonn
Russia	Roszdraznadzor, Address: 4, bld. 1, Slavyanskaya Square, Moscow, 109074
Ukraine	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kiev, 01601
Australia	Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

**Abbreviations:**

CCL – compression class

mmHg – millimetres of a column of mercury (unit for measuring blood pressure in the legs)

**Manufacturer:**

ARIES, a.s., Studenec 309, 512 33 Studenec  
Last reviewed: 21.06.2024

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