## (GB) Instructions for Use

Avicenum PHLEBO 140 FINE knee-high socks  
Avicenum PHLEBO 140 STRONG knee-high socks  
Avicenum PHLEBO 140 FINE thigh-high Stockings  
Avicenum PHLEBO 140 STRONG thigh-high Stockings  
Avicenum PHLEBO 140 FINE Tights

The medical compression stockings that you hold in your hands are a useful device in the treatment of circulatory and lymphatic disorders. You can maximize the benefits of this product by its proper use and maintenance.

**Avicenum PHLEBO 140** stockings are a compression class 1 medical device (I, KT, CCL 1). They were kin and do not contain latex. This compression medical device works by applying precisely defined, made by circular knitting technology. The materials used in our compression medical devices are kind to the sgraduated pressure to the leg, with a gradual reduction in force from the ankle towards the heart. The pressure is greatest in the ankle area at point b (see size chart), reaching values of 18-21 mmHg. A proper fit is necessary to achieve the desired effect. Use the compression stocking when engaged in physical activity; compression products should be removed for longer rest periods (unless otherwise indicated by your doctor). **Avicenum PHLEBO 140** is a medical device intended for use both at home and in medical care facilities.

**Product name:**  
**Avicenum PHLEBO 140 FINE Knee-high socks, product variants:** Closed toe, I,KT, size S normal – L long  
**Avicenum PHLEBO 140 FINE Knee-high socks, product variants:** Open toe, I,KT, size S normal – L long  
**Avicenum PHLEBO 140 STRONG Knee-high socks, product variants:** Closed toe, I,KT, size S normal – L long  
**Avicenum PHLEBO 140 STRONG Knee-high socks, product variants:** Open toe, I,KT, size S normal – L long  
**Avicenum PHLEBO 140 FINE Thigh-high stockings, product variants:** Lace, closed toe, I,KT, size S normal – L long  
**Avicenum PHLEBO 140 FINE Thigh-high stockings, product variants:** Lace, open toe, I,KT, size S normal – L long  
**Avicenum PHLEBO 140 STRONG Thigh-high stockings, product variants:** Hem, closed toe, I,KT, size S normal – L long  
**Avicenum PHLEBO 140 STRONG Thigh-high stockings, product variants:** Hem, open toe, I,KT, size S normal – L long  
**Avicenum PHLEBO 140 FINE Tights, product variants:** Closed toe, I,KT, size S normal – L long

**User profile:**  
**Avicenum PHLEBO 140** is a medical device intended for use both at home and in medical care facilities.

**Intended use of the medical device:**  
This medical device is intended for compression therapy of circulatory disorders and various types of oedema in the leg.

**Indications:** Chronic venous insufficiency – its subjective and objective manifestations, prevention, primary and secondary varicose veins in the legs, superficial thrombophlebitis, oedemas associated with post-thrombotic syndrome after overcoming phlebotrombosis, oedemas and varicose veins in pregnancy, conditions following the occurrence of phlebotrombosis and thrombophlebitis in the legs, Prevention of phlebotrombosis in surgical treatment and subsequently also thromboembolic conditions (especially surgery, orthopaedics, gynaecology), oedemas of the lower limbs (postoperative and post-traumatic, lymphatic, obesity-related venous, idiopathic oedemas of the lower limbs), compression support following sclerotherapy of varicose veins of the legs, compression support following surgical extripation of varicose veins (classic surgery, radiofrequency ablation or laser ablation of varicose veins), travel-related thrombosis (economy class syndrome).

If you are using a compression medical device for the first time, consult your doctor to determine the proper compression class and have regular check-ups during the treatment.

**Contraindications:** acute weeping skin manifestations, acute limb ischaemia, critical limb ischaemia, decompensated cardiac insufficiency, phlegmasia cerulea dolens, phlegmasia alba dolens, hypersensitivity to materials used, limb numbness disorders (e.g. peripheral neuropathy in diabetes mellitus).

Your doctor should consider the use of the compression stockings in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

**Target group:**  
Patients with primary varicose veins, patients after venous thrombosis and patients with the above-mentioned indications. Compression medical stockings are intended for use on healthy, intact skin. This compression medical device is suitable for all age groups and both sexes. There are no restrictions for use during pregnancy.

**Side effects:**  
Side effects may occur, especially if the compression medical device does not fit properly or if it is worn when at rest. There may be bruising, impaired blood flow due to pressure on blood vessels, or impaired sensitivity due to excessive nerve compression. Sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) in places where the compression stocking comes into contact with the skin.

**Fitting your compression stockings:**  
A proper fit is necessary to achieve the desired effect. Use the size chart to determine the correct size by measuring limb circumference at the points indicated.

**Note:**  
Any oedema of the lower limbs should be examined by your doctor. If this is the first time you experience oedema, consult your doctor before using the product. Do not use the product if you have not determined the cause of the oedema. Always consult your doctor before using this product in combination with other medical devices (such as bandages or orthoses). Use a tape measure, take measurements while standing, preferably first thing in the morning – legs may swell over the course of the day, ask someone to help you, measurements may change over the course of treatment – always check your size before using a new pack of compression stockings.

**If the measured values do not correspond to any of the sizes in the table (e.g. point b falls under size S, but point d is size L), visit a specialised medical supplies store or consult the manufacturer. You may need a custom-made compression stocking.**  
**If the measured values differ significantly for the right and left limb, visit a specialised medical supplies store or consult the manufacturer. You may need a custom-made compression stocking.**

**If you are on the border between two sizes, choose the bigger one.** If you are unsure of your size, visit a specialised medical supplies store or consult the manufacturer.  
**If you have difficulty putting on your compression medical stockings, visit a specialised medical supplies store and request a demonstration of the correct method of putting it on.** You can also use one of the aids for putting on the product.  
**If the hem or lace does not stay in place on the leg** (for thigh-high stockings), wipe the silicone strip or studs with an alcohol solution, or shave the hair on the limb at the point of contact of the silicone with the skin.

**With the correctly selected size:** you should feel noticeable but not unpleasant pressure over the entire length of the compression medical device, the compression medical device should not slide down when walking and should not bunch up under the knee or at the instep, the hem or lace edge of the compression stockings should not be too tight (for thigh-high stockings).

Information on compression class, size, shelf life and manufacturer can be found on the label on the medical device packaging.

### Putting on your compression medical stockings:

Pull the compression stocking onto the leg gradually; do not tug on the upper edge alone. Compression medical stockings should be put on in the morning before any physical activity and before any oedema occurs. When pulling on the compression medical stockings, keep them away from sharp objects (e.g. jewellery) or long nails. Increased care of the soles of your feet (rough skin on the heels) will reduce the risk of damaging the product. The compression stockings must be stretched evenly, without folding, to ensure effective pressure distribution. If you use any body cosmetics, allow 15 minutes after application before putting on the medical device. In the event of visible mechanical damage to the product, stop wearing it and replace it with a new one.

### Care and maintenance:

**Avicenum PHLEBO 140** compression medical stockings are intended for repeated use by a single individual. This product may be used safely and effectively for 6 months from the first use if the following guidelines are respected: this compression medical device can be machine washed (max. temperature 30°C, gentle cycle) with a soapy solution or a special detergent for compression medical devices, in a protective cover (washing bag), do not use fabric softener, do not spin-dry, press the wet compression medical device dry between two towels lay flat to dry, away from a direct heat source (radiator, direct sunlight), do not iron compression medical devices, maximum number of washing cycles is 120, store in a dark, dry place in the original packaging, do not alter the compression medical device in any way, do not repair areas damaged by use, do not trim the lace or hem, restore the adhesive quality of the silicone strips on the lace or the silicone nubs on the hem by wiping with an alcohol solution, keep the compression medical device away from organic solvents and bleaching agents

This medical device may be used safely for 60 months from the date of manufacture. This date is part of the product batch number expressed in the format YY/MM-1234567, where YY is the year and MM the month of manufacture.

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