Avicenum **FINE 360 PHLEBO**

(GB) Instructions for Use

Avicenum PHLEBO 360 FINE calf-length Stockings
Avicenum PHLEBO 360 FINE thigh-high Stockings
Avicenum PHLEBO 360 FINE Tights

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Avicenum PHLEBO 360 FINE Wadenstrümpfe
Avicenum PHLEBO 360 FINE Schenkelstrümpfe
Avicenum PHLEBO 360 FINE Strumpfhose

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(GB)

INSTRUCTIONS FOR USE

Dear customer, The medical compression stockings that you hold in your hands are a useful device in the treatment of circulatory and lymphatic disorders. You can maximise the benefits of this product by its proper use and maintenance.

Avicenum PHLEBO 360 FINE stockings are a Compression Class II compression stockings (CCL II). They were made by circular knitting technology. The knitwear structure is lightweight. The materials used in our compression medical stockings are kind to the skin and do not contain latex. This compression medical device works by applying precisely defined, graduated pressure to the leg, with a gradual reduction in force from the ankle towards the heart. The pressure is greatest in the ankle area at point b (see size chart), reaching values of 23–32 mmHg. A proper fit is necessary to achieve the desired effect. Use the stocking when engaged in physical activity; products should be removed for longer rest periods (unless otherwise indicated by your doctor). **Avicenum PHLEBO 360 FINE** is a medical device intended for use both at home and in medical care facilities.

Product name:

Avicenum PHLEBO 360 FINE calf-length Stockings

Product variants:

- open toe, CCL II, size S normal – XXL long
- closed toe, CCL II, size S normal – XXL long

Avicenum PHLEBO 360 FINE thigh-high Stockings

Product variants:

- lace, open toe, CCL II, size S normal – XXL long
- lace, closed toe, CCL II, size S normal – XXL long
- hem, open toe, CCL II, size S normal – XXL long
- hem, closed toe, CCL II, size S normal – XXL long

Avicenum PHLEBO 360 FINE Tights

Product variants:

- open toe, CCL II, size S normal – XXL long
- closed toe, CCL II, size S normal – XXL long

Related accessories:

Each package of **Avicenum PHLEBO 360 FINE** medical stockings, including all variants, contains a textile dressing aid for Avicenum compression stockings. The dressing aid is not a medical device and serves exclusively to assist you in putting on your Avicenum compression stockings. The use of the dressing aid is described below, along with instructions for how best to put on compression stockings.

User profile:

Avicenum PHLEBO 360 FINE is a medical device intended for use both at home and in medical care facilities.

Intended use of the medical device:

This medical device is intended for compression therapy of circulatory disorders and various types of oedema in the leg.

Indications:

- Chronic venous insufficiency – its subjective and objective manifestations, prevention
- Primary and secondary varicose veins in the legs
- Phlebothrombosis
- Superficial thrombophlebitis
- Oedemas associated with post-thrombotic syndrome after overcoming phlebothrombosis
- Oedemas and varicose veins in pregnancy
- Conditions following the occurrence of phlebothrombosis and thrombophlebitis in the legs
- Prevention of phlebothrombosis in surgical treatment and subsequently also thromboembolic conditions (especially surgery, orthopaedics, gynaecology)
- Oedemas of the lower limbs (postoperative and post-traumatic, lymphatic, obesity-related oedemas, idiopathic oedemas of the lower limbs)
- Compression support following sclerotherapy of varicose veins of the legs
- Compression support following surgical extirpation of varicose veins (classic surgery, radiofrequency ablation or laser ablation of varicose veins)
- Venous compression syndromes (Cockett's syndrome – oppression of the left iliac vein by the right-sided iliac artery, oppression of the popliteal vein, oppression of the vein from the surroundings – e.g. by tumour, cyst or surrounding haematoma)
- Venous malformations (Klippel-Trenaunay syndrome, Parkes-Weber syndrome, Maffucci syndrome)
- Lymphedema (primary and secondary)
- Lipoedema
- Therapy of varicose ulcer
- Travel-related thrombosis (economy class syndrome)

If you are using compression medical devices for the first time, it is advisable to consult their use with the attending doctor to determine the compression class or to have the condition checked regularly by a doctor.

Contraindications:

- Acute weeping skin manifestations
- Acute limb ischaemia
- Critical limb ischaemia
- Decompensated cardiac insufficiency
- Phlegmasia cerulea dolens, phlegmasia alba dolens
- Hypersensitivity to materials used
- Limb numbness disorders (e.g. peripheral neuropathy in diabetes mellitus)

Your doctor should consider the use of the compression stockings in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

Target group:

Patients with primary varicose veins, patients after venous thrombosis and patients with the above-mentioned indications. Compression medical stockings are intended for use on healthy, intact skin. Compression medical stockings are suitable for all age groups and both sexes. There are no restrictions for use during pregnancy.

Side effects:

Side effects may occur, especially if the compression medical stockings do not fit properly or they are worn when at rest. There may be bruising, impaired blood flow due to pressure on blood vessels, or impaired sensitivity due to excessive nerve compression. Sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) in places where the compression medical stocking comes into contact with the skin.

How to choose the correct stockings size:

A proper fit is necessary to achieve the desired effect. Use the size chart to determine the correct size by measuring limb circumference at the points indicated.

Note:

Any oedema of the lower limbs should be examined by your doctor. If this is the first time you experience oedema, consult your doctor before using the product. Do not use the product if you have not determined the cause of the oedema. Always consult your doctor before using this product in combination with other medical devices (such as bandages or orthoses).

- use a tape measure
- Take measurements while standing, preferably first thing in the morning – legs may swell over the course of the day
- Ask someone to help you
- Measurements may change over the course of treatment – always check your size before using a new pack of compression stockings


SIZE TABLE:

SIZE	a	b	b1	c	d	e	f	g	h	t
S	19 – 25	20 – 23	25 – 30	32 – 38	31 – 37	33 – 39	41 – 49	47 – 57	110	80
M	22 – 28	23 – 26	28 – 33	35 – 41	34 – 40	36 – 42	45 – 53	52 – 62	120	90
L	25 – 31	26 – 29	31 – 36	38 – 44	37 – 43	39 – 45	49 – 57	57 – 67	120	90
XL	28 – 34	29 – 32	34 – 39	41 – 47	40 – 46	42 – 48	53 – 61	62 – 72	130	100
XXL	31 – 37	32 – 35	37 – 42	44 – 50	43 – 49	45 – 51	57 – 65	67 – 77	130	100

LENGHT			
	short	normal	long
AD	-	34 – 38	38 – 43
AG (AT)	56 – 62	62 – 71	71 – 82

- a foot circumference
- b circumference above the ankle
- b1 circumference 10 cm above the ankle
- c calf circumference
- d circumference below the knee
- e circumference above the knee
- f thigh circumference (in the middle of the thigh)
- g thigh circumference (5 cm under the crotch)
- h hip circumference
- t waist circumference

- AD calf-length stockings
- AG thigh-high stockings
- AT thighs

If the measured values do not correspond to any of the sizes listed in the size chart (e.g. point b corresponds to size S, but point d is already size L), visit a specialised medical supplies store or the manufacturer. It may be that you need a compression medical stockings custom made just for you.

If the measured values differ significantly for the right and left limb, visit a specialised medical supplies store or consult the manufacturer. It may be that you need a compression medical stockings custom made just for you.

If you are on the border between two sizes, choose the bigger one. If you are unsure of your size, visit a specialised medical supplies store or consult the manufacturer.

If you have difficulty putting on your compression medical stockings, visit a specialised medical supplies store and request a demonstration of the correct method of putting it on. You can also use one of the aids for putting on the product.

If the hem or lace does not stay in place on the leg (for thigh-high stockings), wipe the silicone strip or nubs with an alcohol solution, or shave the hair on the limb at the point of contact of the silicone with the skin.

With the correctly selected size:

- You should feel noticeable but not unpleasant pressure over the entire length of the compression medical stockings
- The compression medical device should not slide down when walking and should not bunch up under the knee or at the instep
- The hem or lace of the compression stockings should not be too tight (for thigh-high stockings)

Information about compression class, size, and manufacturer can be found on the tag sewn at the upper edge of the compression stocking. Should the tag be removed in the course of use, sizing information is also knit in on the stocking sole.

Application of compression stockings:

Pull the compression stocking onto the leg gradually; do not tug on its upper edge alone. Compression medical stockings should be put on in the morning before any physical activity and before any oedema occurs. For easier application, use the enclosed Avicenum compression stocking textile dressing aid. The use of the textile dressing aid is illustrated on the enclosed card. When pulling on the compression medical stockings, keep them away from sharp objects (e.g. jewellery) or long nails. Increased care of the soles of your feet (rough skin on the heels) will reduce the risk of damaging the product. The compression stockings must be stretched evenly, without folding, to ensure effective pressure distribution. If you use body cosmetics, wait about 15 minutes after application before putting on the product. In the event of visible mechanical damage to the product, stop wearing it and replace it with a new one.

Care and maintenance:

Avicenum PHLEBO 360 FINE compression medical stockings are intended for repeated use by a single individual. This product may be used safely and effectively for 6 months from the first use if the following guidelines are respected:

- These compression medical stockings can be machine washed (max. temperature 30°C, gentle cycle) with a soapy solution or a special detergent for compression medical devices, in a protective cover (washing bag)
- Do not use fabric softener
- Do not spin-dry, press the wet compression medical stockings dry between two towels
- Lay flat to dry, away from a direct heat source (radiator, direct sunlight)
- Do not iron compression medical stockings
- Maximum number of washing cycles is 120
- Store in a dark, dry place in the original packaging
- Do not alter the compression medical stockings in any way, do not repair areas damaged by use, do not trim the lace or hem
- Restore the adhesive quality of the silicone strips on the lace or the silicone nubs on the hem by wiping with an alcohol solution
- Keep the compression medical stockings away from organic solvents and bleaching agents

These medical devices may be used safely for 60 months from the date of manufacture. This date is part of the product batch number expressed in the format YY/MM-1234567, where YY is the year and MM the month of manufacture.

Care symbols:

When not to use:

Stop using the product immediately if any of the following conditions apply:

- There are signs of mechanical damage – replace with new compression medical stockings
- You notice a deterioration in your health that may be associated with the use of the compression medical stockings – consult your doctor in such case

Composition:

Calf-length stockings: 60% PAD Nylon, 40% elastane LYCRA®

High-high stockings: 60% PAD Nylon, 40% elastane LYCRA®, silicone

Stockings: 60% PAD Nylon, 40% elastane LYCRA®

These medical devices do not contain any active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, or tissues or cells of animal origin or derivatives thereof.

Disposal:

Compression medical stockings can be disposed of with normal household waste or in a textile container. Although the product is intended for use on healthy, intact skin, it may become contaminated with body fluids during use. In such a case, dispose of the compression medical stockings as infectious material.

The material used is non-toxic unless ignited. Never incinerate used compression medical stockings.

Any serious adverse event occurring in connection with the product in question should be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located:

Czech Republic	State Institute for Drug Control, Šrobárova 48, 100 41 Prague 10
Slovak Republic	State Institute for Drug Control, Kvetná 1024/11, 821 08 Bratislava
Poland	Office of Registration of Medical Products, Medical Products and Biological Products, Al. Jerozolimskie 181C, 02-222 Warsaw
Germany	Federal Institute for Drugs and Medical Devices, Kurt Georg Kiesinger Allee 3, 53175 Bonn
Russia	Roszdraznadzor, Address: 4, bld. 1, Slavyanskaya Square, Moscow, 109074
Ukraine	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kyiv, 01601
Australia	Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

List of abbreviations:

CCL – compression class

mmHg – millimetres of mercury column (unit of blood pressure measurement on the leg)

Manufacturer:

ARIES, a.s., Studenec 309, 512 33 Studenec, Czech Republic
Last reviewed: 21.06.2024

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