Avicenum **TRAVEL 360** **PHLEBO**

(GB) Instructions for Use
Avicenum PHLEBO 360 TRAVEL calf-length Stockings

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(GB)

INSTRUCTION MANUAL

Dear customer, in your hands you are now holding compressive medical hosiery that will help you in the treatment of venous or lymphatic diseases. You can maximize the benefits of the product by its proper use and maintenance.

Avicenum PHLEBO 360 TRAVEL is compressive medical hosiery of compression class II (KT II, CCL 2). The product was manufactured by the circular knitting technology. The materials used for the manufacture of the product are non-irritating and latex-free. The mechanism behind the effects of the product is based on the action of precisely defined graduated pressure on the leg, which reduces from the ankle up to the heart. The pressure applied in the ankle area (see the size table) is the greatest and reaches 23 - 32 mmHg. In order to achieve the desired effects, you must select the correct size of the product. The product is intended for physical exercise. You should take off the product for extended periods of rest (unless your doctor determines otherwise).

Avicenum PHLEBO 360 TRAVEL products are medical devices intended for use both at home and in medical facilities.

Product name:

Avicenum PHLEBO 360 TRAVEL Calf-Length Stockings

Product variants:

- closed toe, KT II, sizes S normal - XXL long

User profile:

Avicenum PHLEBO 360 TRAVEL is a medical device intended for use both at home and in medical facilities.

Intended use:

The product is intended for the compression therapy of venous conditions and swelling of the legs of various origin.

Indications:

- chronic venous insufficiency - subjective and objective manifestations, prevention
- primary and secondary varices of the legs
- deep venous thrombosis
- superficial venous thrombophlebitis
- swelling associated with post-thrombotic syndrome after deep venous phlebothrombosis
- swelling and varices during pregnancy
- conditions after phlebothrombosis and thrombophlebitis of the legs
- prevention of deep venous phlebothrombosis in surgery and associated thromboembolic diseases (especially surgery, orthopaedics, gynaecology)
- swelling of the legs (post-operative and post-traumatic, lymphatic, overweight swelling, idiopathic swelling of the legs)
- compression support after sclerotherapy of leg varices
- support after surgical extirpation of varices (standard surgery, radiofrequency ablation, or laser ablation)

- venous compression syndromes (Cockett's syndrome - compression of the left iliac vein by the right iliac artery, compression of the popliteal vein, surrounding compression of veins, e.g., by tumour, cyst, or surrounding haematoma)
- venous malformations (Klippel-Trenaunay syndrome, Parkes-Weber syndrome, Maffucci syndrome)
- lymphedema (primary and secondary)
- lipoedema
- leg ulcer therapy
- travel thrombosis (economy class syndrome)

If this is the first time you use compression medical hosiery, consult your doctor to select the proper compression class of the product or have your medical condition regularly examined.

Contraindications:

- acute weeping skin manifestations
- acute limb ischemia
- critical limb ischemia
- decompensated cardiac insufficiency
- phlegmasia coerulea dolens, phlegmasia alba dolens
- hypersensitivity to materials used
- limb sensitivity disorders (e.g., peripheral neuropathy in diabetes mellitus)

Your doctor should consider the use of the product in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

Target group:

Patients with primary varices, patients after venous thrombosis, and patients with the above indications. The product may be used only on healthy and intact skin. The product is intended for all age groups and both sexes. The product may be used in pregnancy without limitation.

Side effects:

Side effects may occur especially when using a wrong size of the product or if the product is used outside strenuous activities. This may cause pressure sores, impaired blood circulation due to compression of blood vessels, or deterioration of sensitivity due to excessive compression of the nerves. Sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) at the place of contact of the product materials with the skin.

How to choose the correct size:

In order to achieve the desired effects, you must select the correct size of the product. Determine the size according to the size table by measuring the circumference at the points indicated in the figure.

Note:

Any leg swelling should be examined by your doctor. If this is the first time you experience the swelling, consult your doctor before using the product. Do not use

the product if you have not determined the cause of the swelling. Always consult your doctor before combining the product with other medical devices (e.g., sleeves or braces).

- use a tape measure
- take the measurements standing up, preferably in the morning right after waking up - the legs may swell during the day
- get assistance from another person
- measured values may change during the therapy - always verify the correct size before using a new product


SIZE TABLE:

SIZE	a	b	b1	c	d
S	19 – 25	20 – 23	25 – 30	32 – 38	31 – 37
M	22 – 28	23 – 26	28 – 33	35 – 41	34 – 40
L	25 – 31	26 – 29	31 – 36	38 – 44	37 – 43
XL	28 – 34	29 – 32	34 – 39	41 – 47	40 – 46
XXL	31 – 37	32 – 35	37 – 42	44 – 50	43 – 49

LENGTH		
	normal	long
AD	34 – 38	38 – 43

- a foot circumference
 b circumference above the ankle
 b1 circumference 10 cm above the ankle
 c calf circumference
 d circumference below the knee

AD calf-length stockings

If the measured values do not correspond to any of the sizes in the table (e.g., point b corresponds to size S but point d is already size L), visit a specialized medical supplies store or the manufacturer. It could be that you need a custom-made product.

If the measured values differ significantly for the right and the left leg, visit a specialized medical supplies store or directly the manufacturer. It could be that you need a custom-made product.

If you are on the border between two sizes, choose the bigger one. If you have difficulties choosing the size of the product, visit a specialized medical supplies store or directly the manufacturer.

If you have difficulties putting on the product, visit a specialized medical supplies store and request a demonstration of the correct method of applying it. You can also use one of the aids for putting on the product.

With a correctly selected size:

- you should feel noticeable but not unpleasant pressure over the entire length of the product
- the product does not slide nor roll down

Information on the compression class, size, and manufacturer can be found on the tag sewn at the top edge of product. If the tag is removed over the course of using the product, its size is still indicated on the sole of the product.

Method of use:

Put on the product by gradually pulling it over the leg (do not pull only the top edge). It is preferable to put on the product in the morning before any physical strain or formation of swelling. When putting on, the product may not come into contact with sharp objects (e.g., jewellery) or long nails. Intensive care for your feet (no rough skin on the heels) will decrease the risk of damaging the product and increase its life. Stretch the product evenly without creases for effective distribution of pressure. If you use body cosmetics, wait about 15 minutes after application before you put on the product. In the case of visible mechanical damage to the product, stop wearing it and preferably replace it with a new product.

Care and maintenance:

Avicenum PHLEBO 360 TRAVEL is intended for repeated use by a single person. The product will maintain its properties and safety for six months from the first use under these conditions:

- the product can be washed in a washing machine (maximum washing temperature 30 °C, gentle cycle) in a soapy solution or using a special preparation for washing compression hosiery using a protective container (wash bag)
- do not use fabric softeners
- do not use centrifugation; squeeze the wet product between two towels
- dry the product in a horizontal position away from a direct heat source (radiator, sun)

- do not iron the product
- the maximum number of washing cycles is 120
- store in a dry and dark place, preferably in the original packaging
- do not interfere mechanically in the product, do not repair places damaged by use, and do not cut
- the product may not come into contact with organic solvents or bleaching agents

The product is safe to use for 60 months from the date of manufacture. This date is included in the product's batch number in the format YY/MM-1234567 where YY is the year and MM the month of manufacture.

Care symbols:

Prohibited use:

Stop using the product immediately if:

- if it shows signs of mechanical damage (replace the product)
- your health has deteriorated which could be associated with the use of the product - consult your doctor

Composition:

82 % PAD Nylon, 18 % elastan LYCRA®

This medical device does not contain an active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, tissues or cells of animal origin or derivatives thereof.

Disposal:

The product can be disposed of with normal household waste or in a textile container. Although the product is intended for use on healthy and intact skin, it may be contaminated with body fluids during use. In that case, dispose of the product as an infectious material.

The material used is not toxic unless burned. Never dispose of the product by burning it.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Abbreviations:

KT, CCL – compression class

mmHg – millimetres of a column of mercury (unit for measuring blood pressure in the legs)

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